



fitbag[®]
A REVOLUTIONARY MEDIA BAG



B I O S I D E



COMPACT SIZE

LONG SHELF LIFE 24 MONTHS
STERILIZATION VIA GAMMA
RAYS

ZERO RISK OF CONTAMINATION




SPACE SAVING



5 Average Mediabag - 22,5 Lt.

50 fitbag® - 225 Lt.

EASY TO USE



ONCE FULL YOU CAN USE IT **WHENEVER** YOU WANT



Average workflow - 35 minutes



fitbag[®] workflow - 4 minutes

RAPID
RESTORE AND
DISSOLVEMENT

ENERGY SAVING

Transport and storage at room temperature. Empty bags get disposed of with usual laboratory waste.

COVENIENT AND SAFE

No weighing of powder or supplements, no handling of glassware. No need for expensive equipment and maintenance.

QUALITY CERTIFIED

Zero risk of operator error, no risk of overheating (improper sterilization). Compliant with ISO 11133 and ISO 7218.

EVERY USE

fitbag®	1 Litre	2.7 Litres	4.5 Litres	On demand
BUFFERED PEPTONE WATER (BPW)		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
FRASER BROTH	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>
HALF FRASER BROTH		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
MAXIMUM RECOVERY DILUENT (MRD)		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
TRYPTIC SOY BROTH (TSB)	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>
TETRATHIONATE BROTH (TT)	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>
UNICUM BROTH (UB)		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
OTHERS MEDIA	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
PACKAGING	100 x 1 Litres	30 x 2.7 Litres	20 x 4.5 Litres	On demand

B I O S I D E

Via Albert Einstein, snc
Località Cascina Codazza
c/o Parco Tecnologico Padano
26900 Lodi (LO) ITALY

☎ +39 0371 438 498
☎ +39 0371 438 498
✉ info@bioside.it

www.bioside.it